



















T-Cup[®] Compact Multi-Drug Urine Test Cup





















- 2-20 configurations available
 - Cost-effective
 - Smaller, compact design
 - Low shipping costs
 - CLIA-waived
- Same technical features as T-Cup[®]



CLIA Waived Drugs & Cutoff Levels

	Amphetamine (AMP)	1000ng/mL 500ng/mL		Tricyclic Antidepressants (TCA)	1000ng/mL
	Barbiturates (BAR)	300ng/mL		Methamphetamine (MET/mAMP)	1000ng/mL 500ng/mL
	Benzodiazepines (BZO)	300ng/mL		Methylenedioxymethamphetamine (MDMA)	500ng/mL
	Buprenorphine (BUP)	10ng/mL		Morphine (MOP/OPI300)	300ng/mL
	Cocaine (COC)	300ng/mL 150ng/mL		Opiate (OPI)	2000ng/mL
	Cannabinoids (THC)	50ng/mL		Oxycodone (OXY)	100ng/mL
	2-ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP)	300ng/mL		Phencyclidine (PCP)	25ng/mL
	Methadone (MTD)	300ng/mL		Propoxyphene (PPX)	300ng/mL

Forensic Use Panels

	Monoacetylmorphine (6-MAM)	10ng/mL		Synthetic Cannabinoid (K2)	50ng/mL
	Amphetamine (AMP)	300ng/mL 500ng/mL		Ketamine (KET)	300ng/mL 1000ng/mL
	Benzodiazepines (BZO)	100ng/mL 200ng/mL		Kratom (KRA)	300ng/mL
	Cocaine (COC)	100ng/mL 150ng/mL		Methamphetamine (MET or mAMP)	300ng/mL 500ng/mL
	Cotinine (COT)	200ng/mL		Morphine (MOP or OPI)	100ng/mL
	EDDP	100ng/mL		Methadone (MTD)	200ng/mL
	Ethyl Glucuronide (EtG)	500ng/mL		Methaqualone (MQL)	300ng/mL
	Fentanyl (FTY)	20ng/mL		Cannabinoids (THC)	25ng/mL 40ng/mL
	Hydromorphone (HMO)	300ng/mL		Tramadol (TRA)	100ng/mL 200ng/mL 1000ng/mL

Generic stock adulteration testings include Creatinine (CR), Specific Gravity (S.G.) and pH. Other adulterates like Glutaraldehyde (GL), Oxidant (OX) & Nitrite (NI) are available through special order.

※ **Kit Size / Case Size**
25 / 4

Item Code	Description	Complexity
BB-UD01	12-Panel T-Cup Compact: AMP, BAR, BUP, BZO, COC, mAMP, MDMA, MOP, MTD, OXY, PCP, THC	CLIA WAIVED