

AdvanCare™ Nitrile Exam Gloves

Features

- Advanced formulation for exceptional comfort and fit providing better tactile sensitivity
- 100-count & 200-count packaging for various choices
- Tested for use with over 60 chemicals and approved for use with chemotherapy drugs*




Quality Standards

- Complies with ASTM 6978



Applications

-  Medical Purpose / Examination
-  Healthcare and nursing
-  Industrial purpose / PPE
-  General housekeeping
-  Laboratory
-  IT Industry



SKU	Size	Color	Package	Box Size	Carton Dimension
ANBM1001*	XS-XL	Blue	100 pcs/box, 10 boxes/ctn	230*125*60mm	315*258*245mm
ANBM2001*	XS-XL	Blue	200 pcs/box, 10 boxes/ctn	245*130*80mm	415*270*258mm
ANBM1003*	XS-XL	Violet Blue	100 pcs/box, 10 boxes/ctn	230*125*60mm	315*258*245mm
ANBM2003*	XS-XL	Violet Blue	200 pcs/box, 10 boxes/ctn	245*130*80mm	415*270*258mm



* The last "*" means Size: 3-XS 4-S 5-M 6-L 7-XL.



BigBen Medical Supplies
 1420 Washington Boulevard Detroit, MI
 Call center No: +1 3135664443
 www.Bigbenmedical.com