


# BigBen Medical Supplies



## TECHNICAL DATA SHEET

Type	Disposable Nitrile Gloves, Powder Free, Non Sterile, Ambidextrous, 9"														
Grade	Medical Grade			Surface	Finger Textured			Glove Length(mm)	>=230 ASTM D6319						
Material	NBR			Packing	100pcs/box;10boxes/ctn			Weight(g)	M=3.5g+/-0.3						
Cuff	Rolled-Beaded			Shelf Life	5 Years			Powder Level	≤2mg/glove						
Product	Color	Size	Weight(g)	Physical Dimension					Physical Property						
				Length (mm)	Palm Width (mm)	Cuff Thickness (mm)	Palm Thickness (mm)	Finger Thickness (mm)	Force at Break(Mpa)		Elongation		Water Leak AQL	Appearance AQL	
Disposable Nitrile Gloves		S	3.0+/-0.3	≥230	80±10	0.07±0.03	0.08±0.03	0.110±0.03	≥14	≥14	≥500	≥400	1.5	Critical AQL1.0 Major AQL2.5 Minor 4.0	
		M	3.5+/-0.3	≥230	95±10	0.07±0.03	0.08±0.03	0.110±0.03	≥14	≥14	≥500	≥400	1.5	Critical AQL1.0 Major AQL2.5 Minor 4.0	
		L	3.8+/-0.3	≥230	110±10	0.07±0.03	0.08±0.03	0.110±0.03	≥14	≥14	≥500	≥400	1.5	Critical AQL1.0 Major AQL2.5 Minor 4.0	
		XL	4.1+/-0.3	≥230	120±10	0.07±0.03	0.08±0.03	0.110±0.03	≥14	≥14	≥500	≥400	1.5	Critical AQL1.0 Major AQL2.5 Minor 4.0	
No.	sort		Critical defects					Major defects			Minor defects				
1	torn	big tearing	√					-			-				
		scratch	≥10mm					<10mm			-				
		visible hole	≥10mm					<10mm			-				
2	dirt	greasy dirt	≥10mm,serious					5-10mm,visible			<5mm,slightly				
		yellow dirt	√					-			-				
3	black/PU spots,impurities		uncountable					≥3mm or 0.5-3mm,>3/side			0.5-3mm,≤3/side				
4	beading		broken;torn;none-beading					impurity;knot;black;bright PU ≥ 3mm;inward ≥ 10mm; minuteness(≤0.5cm)			bright PU 0.5-3mm; inward<10mm;thin 0.5-0.8				
5	sticky		-					sticky cuff severely			slightly				
6	Extreme color variation		entirely,severely					finger,cuff,severely			cuff,slightly,unconspicuous				
7	Bad smell		-					severely irritating smell			slight odor				
8	Embedded foreign particles		insect,hair,etc.					-			-				