

## Urinalysis Tests are now available Conveniently Dropship from Houston

Utilize a series of visual, semi-quantitative tests simultaneously to analyze urine properties that aid in detecting/monitoring conditions such as:

- Diabetes
- Metabolic Abnormalities
- Liver Disease
- Kidney Function
- Urinary Tract Infection

Results read within 2 minutes. These tests can be used along with other diagnostic information to rule out certain disease states and determine if further microscropic analysis is needed.

## Testing Parameters

Urobilinogen Protein Glucose
Bilirubin Nitrite Specific Gravity
Ketone Ascorbic Acid pH
Blood Leukocytes





510(k) cleared, CLIA Waived Prescription Use Only

