

**When You Can't
Afford to Wait...**

Status Stik™

Instant Test for Drugs of Abuse.


IMMEDIATE
results in 5 minutes.

SIMPLE
*one step dip is easily and accurately
performed by anyone.*

CONVENIENT
reduced handling of urine.

COMPREHENSIVE
*test menu allows you to choose
which drugs to test.*

97-99% CORRELATION
*to GC/MS, FDA cleared, and standardized
to NIDA (SAMHSA) cutoff levels.*



Status Stik™

Accurate results in minutes to enable timely, effective decisions.

RAPID DRUGS OF ABUSE TEST

- Put new hires to work immediately.
- Get timely results in probable cause and post-accident situations.
- Protect your employees from drug abuse.
- Reduce employee downtime while waiting for drug test results.
- Increase productivity
- Lower drug testing costs associated with outside labs.

Procedure

Remove cap and immerse until specimen starts migrating through result window. Read results in 5 minutes.*



Negative -

a visible line for drug tested

Positive -

no visible line for drug tested

*see package insert for details



Reading the Test

Ordering Information

| Description | 10 Test Kit Stock No. | 35 Test Kit Stock No. |
|---------------------------|-----------------------|-----------------------|
| THC/COC Panel | 12110 | 12135 |
| THC/COC/OPI Panel | 13010 | 13035 |
| THC/COC/MET Panel | 13110 | 13135 |
| THC/OPI/COC/AMP Panel | 14010 | 14035 |
| THC/OPI/COC/MET Panel | 14110 | 14135 |
| THC/OPI/COC/AMP/PCP Panel | 15010 | 15035 |
| THC/OPI/COC/MET/PCP Panel | 15110 | 15135 |

Free on-site training and assistance are available from LifeSign's professional sales staff and technical service departments. For more information or to contact your local sales representative call LifeSign at

1-800-526-2125 or visit us on line...

@ www.lifesignmed.com