

# Status hCG



Rapid immunassay for the qualitative detection of human chronic gonadotropin in serum and urine.

## Status hCG Cassette

For the Detection of hCG in Serum or Urine.

- **CLIA WAIVED** for urine
- Serum and urine 25mIU/ml sensitivity
- Results in 3-5 minutes
- Minimal susceptibility to the Hook effect by hCG  $\beta$  core fragment



## Status hCG Strip

For the Detection of hCG in Urine.

- **CLIA WAIVED**
- Simple "dip and read" procedure
- Results in 3-5 minutes
- Minimal susceptibility to the Hook effect by hCG  $\beta$  core fragment

	CPT CODE
Urine	81025QW
Serum	84703

ORDERING INFORMATION		
PRODUCT	# OF TESTS	ITEM #
Status hCG Strip	35	35235
Status hCG Urine	35	35135
Status hCG Combo Urine/Serum	35	35035
Status hCG Urine Control	2 x 2.5 ml	135-0000A

Proudly made in the USA

For more information or to place an order, contact:



A PBM Group Company  
 85 Orchard Road • Skillman • New Jersey • 08558  
 800.526.2125 • 732.246.3366 • Fax: 732.246.0570 • www.lifesignmed.com

Printed in USA.  
 All rights reserved.  
 78908 - 05/01/2019