

HIV Product Catalogue



Ordering Information

| Product | Catalog # |
|-----------------------------------------------------------|-----------|
| SURE CHECK® HIV 1/2 Assay | 60-9507-0 |
| STAT-PAK® HIV 1/2 Assay | 60-9505-1 |
| DPP® HIV 1/2 Assay | 65-9500-0 |
| DPP® HIV 1/2 Rapid Test Control Pack | 60-9552-0 |
| HIV Rapid Test Control Pack (for STAT-PAK® & SURE CHECK®) | 60-9549-0 |
| HIV Control Accessory Kit | 60-9549-1 |

Reimbursement Codes: 86703 or 86703QW
G0435 (for Medicare)

Performance*

| Product | CLIA WAIVED | | | | |
|-------------------------------------|-------------|-------------------|--------------------|--------|-------|
| | Oral Fluid | Fingerstick Blood | Venous Whole Blood | Plasma | Serum |
| SURE CHECK® HIV 1/2 FDA Approved | Sensitivity | N/A | 99.7% | 99.7% | 99.7% |
| | Specificity | N/A | 99.9% | 99.9% | 99.9% |
| STAT-PAK® HIV 1/2 FDA Approved | Sensitivity | N/A | 99.7% | 99.7% | 99.7% |
| | Specificity | N/A | 99.9% | 99.9% | 99.9% |
| DPP® HIV 1/2 FDA Approved | Sensitivity | 98.9% | 99.8% | 99.9% | 99.9% |
| | Specificity | 99.9% | 100% | 99.9% | 99.9% |

* Refer to product insert for details.



Visit: www.chembio.com
Call: 1.844.CHEMBIO (844.243.6246)
Email: info@chembio.com

© Copyright 2019, Chembio Diagnostic Systems, Inc. The Chembio logo, SURE CHECK, STAT-PAK, and DPP are trademarks of Chembio Diagnostic Systems, Inc.
® Indicates a Chembio trademark registered in the USA. MS-19-008 Rev 1

Rapid HIV Testing Solutions

SURE CHECK® HIV 1/2 Assay



World's Smallest Sample Volume

- 2.5 µl

Easy to Use

- All-in-one sample collection and test device

Convenient

- Each test includes safety lancet and bandage



SIMPLE

STAT-PAK® HIV 1/2 Assay



Small Sample Volume

- 5 µl

Superior Performance

- Only test that meets the WHO recommended thresholds for performance*


Convenient

- Also CLIA waived for venous whole blood



RELIABLE

DPP® HIV 1/2 Assay



Low Sample Volume


- 10 µl

Flexible

- Works with oral fluid, blood, serum, plasma

Convenient

- Safe closed system for sample handling



FLEXIBLE