

HIV Product Catalogue



Ordering Information

Product	Catalog #
SURE CHECK® HIV 1/2 Assay	60-9507-0
STAT-PAK® HIV 1/2 Assay	60-9505-1
DPP® HIV 1/2 Assay	65-9500-0
DPP® HIV 1/2 Rapid Test Control Pack	60-9552-0
HIV Rapid Test Control Pack (for STAT-PAK® & SURE CHECK®)	60-9549-0
HIV Control Accessory Kit	60-9549-1

Reimbursement Codes: 86703 or 86703QW
G0435 (for Medicare)

Performance*

Product	CLIA WAIVED				
	Oral Fluid	Fingerstick Blood	Venous Whole Blood	Plasma	Serum
SURE CHECK® HIV 1/2 FDA Approved	Sensitivity	N/A	99.7%	99.7%	99.7%
	Specificity	N/A	99.9%	99.9%	99.9%
STAT-PAK® HIV 1/2 FDA Approved	Sensitivity	N/A	99.7%	99.7%	99.7%
	Specificity	N/A	99.9%	99.9%	99.9%
DPP® HIV 1/2 FDA Approved	Sensitivity	98.9%	99.8%	99.9%	99.9%
	Specificity	99.9%	100%	99.9%	99.9%

* Refer to product insert for details.

Rapid HIV Testing Solutions

SURE CHECK® HIV 1/2 Assay



World's Smallest Sample Volume

- 2.5 µl

Easy to Use

- All-in-one sample collection and test device

Convenient

- Each test includes safety lancet and bandage



SIMPLE

STAT-PAK® HIV 1/2 Assay



Small Sample Volume

- 5 µl

Superior Performance

- Only test that meets the WHO recommended thresholds for performance*


Convenient

- Also CLIA waived for venous whole blood



RELIABLE

DPP® HIV 1/2 Assay



Low Sample Volume


- 10 µl

Flexible

- Works with oral fluid, blood, serum, plasma

Convenient

- Safe closed system for sample handling



FLEXIBLE